

The logo on the left consists of a white cross inside a circle, with horizontal lines passing through the cross and the circle. The text 'Medioplast' is in a large, bold, white sans-serif font, and 'Plasters' is in a smaller, bold, white sans-serif font below it. The background is a vibrant green with diagonal yellow and white lines, and a blue, angular, crystalline shape in the lower center.

# Medioplast

## Plasters



The overriding aim of Neomedic is to continue to maximize customer value through an unrelenting push for quality that meets our stated benchmarks and exceeds our customer's expectations.

We must at all times demonstrate a commitment to quality in each of the following fields:

**Commitment to a Quality Product:**

- By understanding our ethical responsibilities as vendors of medical products.
- By selecting and selling products that meet recognised quality benchmarks in each of their fields.
- By accepting that our actions have a direct effect on human lives.

**Commitment to a Quality Service:**

- By understanding that there are certain minimum standards for the quality of communication in terms of tone, content and feedback that all our customers have a right to expect from us.
- By acknowledging and answering all customers communication on the day it is received.
- By ensuring that all orders are delivered on time or by providing valid reasons for failing to do so.

**Commitment to a Competitive Price:**

- By ensuring that our prices are most competitive for products of equivalent quality.
- By re-iterating our commitment to better healthcare through ethical pricing.

By committing ourselves to these quality statements we will achieve our simple aim of building a strong company that adds value to every transaction.





# about Mediplast Plasters

Mediplast is a range of high quality Roll and Strip plasters designed to meet the vast majority of Clinical and First Aid requirements.

The Mediplast range consists of 5 different Roll plaster types in both Premium (snap and spool ring) and Hospital (wax-paper wrapped) packs.

The Strip plaster range combines tough highly adhesive plasters, with gentle and discreet plasters. There is also a range of strip plasters uniquely designed to appeal to children. This range is further complemented by the 1 metre Dressing Strips that allow greater flexibility of wound treatment.

All Mediplast plasters are manufactured according to ISO 9002 and EN 46002 standards and are CE certified. In order to enhance usability, all Mediplast packaging is colour coded for size and through the use of high quality illustrations; user instructions have never been easier to follow.

Mediplast is therefore a comprehensive, high quality range of plasters, full of value added features and presented in exquisite packaging.



# Mediplast Strip Plasters

## Fabric Plasters

Mediplast Fabric plasters are made from tough, high performance materials and are designed to promote the efficient healing of small wounds, cuts, scrapes, and other minor injuries.

- Quality materials stretch and conform to different shapes and forms. Guarantees freedom of movement and activity.
- Larger, softer pad with antiseptic; cushions, protects and aids in faster healing.
- Powerful adhesive seals and prevents contamination of wound whilst keeping plaster firmly in place.
- Fabric Weave ventilates wound and promotes, quick natural healing.

## Waterproof Plasters

Mediplast Waterproof plasters are made from strong yet flexible material that stretch without misshaping and keeps the wound sealed at all times.

- Flexible material stretches with every movement without misshaping.
- Larger, softer pad with antiseptic; cushions, protects and aids in faster healing.
- Special adhesive ensures plaster stays on, even when in water.
- Micro-holes allow air in for natural healing, whilst keeping water and germs out.

## Hypoallergenic Plasters

The finest, smoothest materials and low-allergy adhesives have been combined to create a plaster that is recommended for use on children and sensitive skin areas. Mediplast Hypoallergenic plasters are the caring choice for wound treatment of delicate skins.

- Soft and stretchable materials conform to every movement.
- Larger, softer pad with antiseptic; cushions, protects and aids in faster healing.
- Low-allergy adhesive, ensures that plaster is securely but gently held in place.
- Smooth material weave for sensitive wound protection.

## Clear Plasters

Mediplast Clear Strip plasters allow for the discrete, washproof and efficient healing of small wounds, cuts, scrapes and other minor injuries.

- Clear, discreet material for washproof healing.
- Larger, softer pad with antiseptic; cushions, protects and aids in faster healing.
- Special adhesive ensures plaster stays on even when in water.
- Micro-holes allow air in for natural healing, whilst keeping water and germs out.

## FunStrips Kids Plasters

Mediplast FunStrips Kids plasters are multicoloured waterproof strips for fun based, quick, natural healing and protection of wounds, cuts, scrapes and grazes, even when exposed to the rough and tumble of a child's fun packed day.

- Flexible material stretches with every movement without misshaping.
- Soft pad, cushions and protects wound.
- Special hypoallergenic adhesive ensures plaster stays on for extended periods.
- Micro-holes allow air in for natural healing, whilst keeping water and germs out.

# Mediplast Roll Plaster



## 1 Metre Fabric Dressing Strip

**High quality fabric plaster with special hygienic absorbant antiseptic pad.**

- Quality materials made from premium quality fabric.
- Strong and durable adhesive seals and prevents contamination of wound while keeping plaster firmly in place.
- Special hygienic absorbant pad with antiseptic cushions, protects, and aids in healing.

### USES

- Dressing of large cuts, scrapes and grazes and post incision wound.



## Fabric Plaster

**Traditional flexible fabric plaster with extra-strong adhesive.**

- Quality materials made from premium quality woven cotton.
- Strong, non-stretch and air permeable materials are ideal for securing larger dressings in place while offering a certain amount of flexibility.
- Extra-strong adhesive ensures dressing stays in place over long periods.

### USES

- Use to secure larger dressings in place.





## Zinc Oxide Plaster

**Traditional non-stretch plaster with extra-strong adhesive.**

- Quality materials made from premium grade cotton.
- Strong non-stretch materials are permeable to air and water and are ideal for securing larger dressings in place while offering a certain amount of flexibility.
- Extra-strong zinc oxide adhesive ensures dressing stays in place.

### USES

- Use to secure larger dressings in place.
- Use for strapping joints.
- Use as traction, long term bandage.




## Silk Tape

**High quality silk effect material with hypoallergenic adhesive.**

- High quality silk effect material is ideal for discreet, flexible securing of dressing.
- Bi-directional tearing, easy to tear and write on, radio and x-ray transparent material is ideal for fixing and marking small to medium size dressings.
- Hypoallergenic adhesive for sensitive skin, seals and prevents contamination of wound while keeping plaster firmly in place.

### USES

- Use on regular and bulky fixation of dressings and in sensitive areas.
- Use to secure small tubes, catheters, canulae and other devices.



# Mediplast



## Surgical Paper Tape

**Soft Non-woven Paper Tape. Virtually painless upon removal.**

- Soft Non-Woven paper tape made from hypoallergenic materials.
- Strong yet tearable, radio and x-ray transparent, air and water permeable material prevents collection of moisture and is ideal for securing dressings on delicate and sensitive skin.
- Extra-strong adhesive ensures dressing stays in place.

### USES

- Use to secure small and medium dressings on sensitive and delicate skin.
- Use to secure electrodes and other devices.
- Recommended for use on young and elderly skin.



## Clear First Aid Tape

**High quality clear microperforated waterproof material with hypoallergenic adhesive.**

- Premium quality porous materials allow wounds to breathe and heal efficiently.
- Bi-directional tearing and easy to mark material is ideal for fixing and marking small to medium size dressings.
- Hypoallergenic adhesive for sensitive skin seals and prevents contamination of wound while keeping plaster firmly in place.

### USES

- Use on smaller dressings and in sensitive areas.
- Use to secure small tubes, canulae, catheters and other devices.





 **Neomedic**  
LIMITED

Middlesex  
United Kingdom  
[marketing@neomedic.co.uk](mailto:marketing@neomedic.co.uk)

[www.Neomedic.co.uk](http://www.Neomedic.co.uk)